	Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date:	Workflow Versions
	5/1/2025	PHB v3.0.0

1. PURPOSE/SCOPE

To standardize the process of analyzing SARS-COV-2 (SC2) next generation sequencing (NGS) data using Theiagen's TheiaCoV_ClearLabs_PHB workflow in Terra to generate assemblies, quality control (QC) metrics, and determine Nextclade clade and Pangolin lineage assignments. Acceptable data types include ClearLabs raw read file format. Read the documentation [here](#).

2. REQUIRED RESOURCES

- Computer
- Internet connection: at least 10 and 5Mbps for download and upload speeds, respectively
- Internet browser
 - Google Chrome, Firefox, or Edge
- Google account
- Terra account, linked to Google account
- ClearLabs raw sequencing read files uploaded to Terra workspace, see [TG-TER-03](#)
- Theiagen's TheiaCoV_ClearLabs_PHB workflow in Terra, see [appendix 10.1](#)

IMPORTANT NOTES

- Metadata *column headers* and *workflow input* text indicated in **gray** in this SOP are customizable; **black** is required text
- Terra data table column headers become available as workflow inputs when running workflows, search for them in workflow input dropdowns using the prefix **this** to filter
- Filter for workspace data and files in workflow input dropdowns using the prefix **workspace**.


3. RELATED DOCUMENTS


Document Number	Document Name
TG-TER-03	Getting Started in Terra: Importing Reads, Metadata, Workflows, and More

4. PROCEDURE


4.1 RUNNING THE THEIACOV WORKFLOW

- Open Terra and navigate to the **workflows** tab within the workspace containing SC2 data (Fig 1)
- Select the **TheiaCoV_ClearLabs_PHB** (Fig 1)




	Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date:	Workflow Versions
	5/1/2025	PHB v3.0.0

9. Specify the desired **dataset tags** and **docker image** inputs
 - a. To run TheiaCoV_ClearLabs_PHB for the first time or configure with the newest dataset tags and docker images, upload the TheiaCov input json file on the inputs tab by navigating to the Key Resources Notion page titled Docker Image and Reference Materials for SARS-CoV-2 Genomic Characterization
 - i. **NOTE:** TheiaCoV PHB v2.0.0+ workflows are not backwards compatible with older versions of Nextclade; **use Nextclade Dataset Tag 2024-04-15—15-08-22Z or newer**
 - ii. Expand the **TheiaCoV in PHB (v2.0.0 or higher)** section, followed by the **Terra.Bio Input JSONs for PHB v2.0.0 or higher**; click on the json file associated with the ClearLabs platform, **TheiaCoV_ClearLabs_PHB_2025-04-02.json**, or newer
 - iii. **Right click** and **save** the file (text does not have to be selected to save properly)
 - iv. Return to the workflow in Terra, click **upload json** (Fig 4, red circle), **select** the saved json file, and click **open**
 - b. To run the workflow with previously saved dataset tags and docker images, no changes are needed
 - c. To add docker images and dataset tags as workspace files for availability in input dropdowns, refer to [appendix 10.2](#)



10. Set the first and third attributes in the table to **this.reads** and **this.KK_ClearLabs_Test_id**, respectively (Fig 4)
 - a. Where **this.KK_ClearLabs_Test_id** is the unique name of your data table in Terra
11. Manually choose the **primer_bed** file for the primer set used to sequence samples
 - a. Labs using the Artic V4-1 will choose **workspace.Artic_V4-1_primer_bed**; for other primer bed files, see [Docker Image and Reference Materials for SARS-CoV-2 Genomic Characterization](#) for available primer bed files




**Analyzing SARS-CoV-2 Data in Terra using Theiagen's
TheiaCoV ClearLabs Workflow**

Document TG-SC2-CL, Version 4

Date:	Workflow Versions
5/1/2025	PHB v3.0.0


b. To add workspace files for availability in input dropdowns, refer to [appendix 10.2](#)




12. Specify outputs by clicking on the **outputs** tab and **use defaults** (Fig 5)
13. Click **save**
14. Launch the workflow by clicking **run analysis** (Fig 5); enter desired comments and click **launch**

4.2 QUALITY ASSESSMENT OF THEIACOV OUTPUTS

1. Navigate to the **data** tab of the workspace containing SC2 data and open the pertinent data table
2. Click **Select Columns** (Fig 6, green rectangle) and select **none** to deselect all output columns (Fig 6, yellow highlight)
3. To simplify the table, select the three following outputs that will be used to make a QC assessment: **assembly_length_unambiguous**, **Number_N**, and **percent_reference_coverage**
 - a. *Optional:* save this selection by clicking in the **save this column selection** field and naming it (e.g. **QC_assessment**); *do not include any spaces* in the name (Fig 6, red rectangle)
 - b. Click **done**




	Analyzing SARS-CoV-2 Data in Terra using Theiagen's	
	TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date: 5/1/2025	Workflow Versions PHB v3.0.0

4. Optional: add a column to record QC PASS/FAIL by clicking **edit**, **add a column** (Fig 7)

a. Name the new column (e.g. *QC_Call*); **do not include any spaces**

b. Set the value type as **string**

c. Click **save**




5. Use table 1 to assess the quality of each sample's genome assembly (see next page) &/or lab-specific quality metrics
6. Optional: notate in the *QC_assessment* field for each sample **PASS** or **FAIL** by clicking the pencil icon in the corresponding field (Fig 7, red circle)
7. For samples that pass the guidance thresholds, proceed to [section 4.3](#)
- a. For samples that do not pass guidance thresholds, resequence
- i. Samples not meeting guidance thresholds indicated here may proceed to analysis at the discretion of the laboratory

Table 1. Guidance thresholds for genome assembly QC

QC Metric	Data Table Column	Guidance Threshold ¹
Number N	<i>number_N</i>	<5kbp
Assembly length unambiguous	<i>assembly_length_unambiguous</i>	>24kbp
Percent reference coverage	<i>percent_reference_coverage</i>	>83%

¹ Metrics and thresholds presented for guidance only as there are currently no standard assembly metric requirements; internal validation procedures will ultimately define acceptable assembly QC parameters


	Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date:	Workflow Versions
	5/1/2025	PHB v3.0.0

4.3 DETERMINING SARS-CoV-2 CLADES, LINEAGES, AND WHO VARIANTS OF CONCERN (VoC)

1. Navigate to the **data** tab of the Terra workspace containing SC2 data of interest
2. **Open the data table** by clicking on the name of the data table in the left sidebar
3. View **settings** above the data table (Fig 6), select **none** (Fig 7)
4. Select the following columns: **nextclade_clade** and **pango_lineage**
 - a. *Optional: save this column group for future use by clicking the **save this column selection** field, naming it (e.g. **SC2_Results**), and clicking **save***
5. Click **done**
6. Determine the Nextclade clade for each sample
 - a. In the data table, find the column titled **nextclade_clade**; result formats will use the following nomenclature: **21L (Omicron)** where:
 - i. **21L** indicates the sample clade and
 - ii. In parentheses, **(Omicron)**, contains the WHO variant of concern classification
 1. *Not every sample will belong to a WHO classification*
 - b. *Samples indicated as recombinant may indicate a case where multiple strains have combined during viral replication producing a new lineage*
 - c. *More information on SARS-CoV-2 recombinants can be found at the following Github site: [pipeline-resources/docs/sc2-recombinants.md at main · pha4ge/pipeline-resources · GitHub](https://github.com/pha4ge/pipeline-resources/blob/main/docs/sc2-recombinants.md)*
7. Identify the Pangolin lineage for each sample
 - a. In the data table, find the column titled **pango_lineage**; nomenclature will be similar to the following: B.1.167
 - b. *For more information on each of the lineages, visit https://cov-lineages.org/lineage_list.html*
8. Follow lab-specific QC, resulting, and reporting procedures, as applicable

5. QUALITY RECORDS

- Raw read files
- Workflow version and input parameters
- Reference sequence, if applicable
 - a. SC2: Wu, F., et al. (2020). Severe acute respiratory syndrome coronavirus 2 isolate Wuhan-Hu-1, complete genome. NC_045512.2. [FASTA Genome Assembly]. NCBI. <https://www.ncbi.nlm.nih.gov/nuccore/1798174254>.
- Sample read, assembly, and result-specific QC metrics
- All workflow outputs relevant to results, including tool and database versions

	Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date:	Workflow Versions
	5/1/2025	PHB v3.0.0

6. TROUBLESHOOTING

- Consult with internal staff familiar with this procedure or contact support@theiagen.com for troubleshooting inquiries
- For document edit requests, contact support@theiagen.com

7. LIMITATIONS


- This SOP is written for the analysis of SC2 data; v2+ of the TheiaCoV_ClearLabs_PHB workflow is also compatible with the following pathogens: monkeypox virus (MPXV), human immunodeficiency virus (HIV), west nile virus (WNV), influenza virus, and respiratory syncytial viruses A and B (RSV). Refer to [Theiagen Public Health Resources documentation](#) for organism-specific parameters and details.
- TheiaCoV PHB v2.0.0+ workflows are not backwards compatible with older versions of Nextclade; use Nextclade Dataset Tag **2024-04-15—15-08-22** or newer*

8. REFERENCES

- Smith, E., Wright, S., & Libuit, K. (2022, June 28). *Identifying SARS-CoV-2 Recombinants*. Github. Retrieved June 16, 2023, from <https://github.com/pha4ge/pipeline-resources/blob/main/docs/sc2-recombinants.md#identifying-sars-cov-2-recombinants>
- O'Toole, Áine et al. "Tracking the international spread of SARS-CoV-2 lineages B.1.1.7 and B.1.351/501Y-V2 with grinch." *Wellcome open research* vol. 6 121. 17 Sep. 2021, doi:10.12688/wellcomeopenres.16661.2
- Libuit, Kevin G., Emma L. Doughty, James R. Otieno, Frank Ambrosio, Curtis J. Kapsak, Emily A. Smith, Sage M. Wright, et al. 2023. "Accelerating Bioinformatics Implementation in Public Health." *Microbial Genomics* 9 (7). <https://doi.org/10.1099/mgen.0.001051>
- Theiagen Genomics [Public Health Bioinformatics Workflow Documentation](#)



9. REVISION HISTORY

Revision	Version	Release Date
Document creation	1	7/2023
Added TG-TER-04 reference, uncheck call caching, updated input json, figures, and formatting	2	9/2023
Removed section 4.1 for creating a metadata tsv file (refer to TG-TER-03 and TG-TER-04 for details); updated quality records and limitations sections; added primer bed file upload instructions; added appendices 10.1 and 10.2	3	5/2024
Updates for version release, aligning with Terra interface, formatting	4	5/2025


	Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date:	Workflow Versions
	5/1/2025	PHB v3.0.0


10. APPENDICES


10.1 IMPORTING THE THEIACOV_CLEARLABS_PHB WORKFLOW FROM DOCKSTORE




1. In the **Terra workspace** of interest, open the **workflows** tab and click **find a workflow** (Fig 8)
2. In the pop-up window, click **dockstore** (Fig 9)
3. In the top banner click **Organizations**; then click **Theiagen Genomics** (Fig 10)
4. Open the **Public Health Bioinformatics (PHB)** collection (Fig 11)



	Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date:	Workflow Versions
	5/1/2025	PHB v3.0.0




5. To find the TheiaCoV_ClearLabs_PHB workflow in Windows environments, hold **Ctrl + F** and search **TheiaCoV_ClearLabs**, then click on the link (Fig 12)




6. Click **Terra** to launch the workflow in Terra (Fig 13)



7. Choose the **destination workspace** in the dropdown and click **import** (Fig 14)






Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
Document TG-SC2-CL, Version 4	
Date:	Workflow Versions
5/1/2025	PHB v3.0.0




10.2 ADDING WORKSPACE DATA ELEMENTS

1. Navigate to the Terra workspace where analysis will be run
2. To upload local files, open the Files tab in the right-side panel of the workspace (Fig 15, box)
 - a. Click upload (Fig 16)
 - b. Once the upload is complete, right click on the file name and click copy link



	Analyzing SARS-CoV-2 Data in Terra using Theiagen's TheiaCoV ClearLabs Workflow	
	Document TG-SC2-CL, Version 4	
	Date:	Workflow Versions
	5/1/2025	PHB v3.0.0

- Open the **workspace data** tab (Fig 15, circle)
- Click Edit and Add variable in the top tool bar (Fig 17)



- Click in the **key field** and **name the element** being added (Fig 18)
 - E.g. to add the Artic v4-1 primer bed file, the key **Artic_v4-1_primer_bed** may be used
- In the value field, choose **string** as the value type
 - Paste the file path**; the value should start with **gs://**
 - NOTE:** For other string elements like dataset tags and docker images **paste the ID value**
 i.E.g. for the nextclade docker image, add **nextstrain/nextclade:2.14.0**
 ii. Always ensure the docker images and dataset tags are aligned with versions used for internal validation procedures
- Optional:** A description may be added to denote the date updated with staff initials
- Click the blue check mark on the right-hand side of the variable to save it
 - The variable will now be available as a workflow input which can be found by typing the prefix **workspace.** plus the key name **artic_v4-1_primer_bed**
 i. e.g. **workspace.artic_v4-1_primer_bed**

Key	Value	Description
Artic_V4-1_primer_bed	V4-1_nCoV-2021.primer.bed	
Artic_V4_primer_bed	V4_nCoV-2021.primer.bed	
FreyjaLineageMetadata	curated_lineages.json	Taken from Freyja_Workflows Demo Data
FreyjaUsherBarcodes	usher_barcodes.csv	Updated 8/3/23; taken from Freyja_Workflows Demo Data
Freyja_ReferenceGenome	nCoV-2019.reference.fasta	MN908947.3
Midnight_primer_bed	Midnight_PrimerSARS-CoV-2.scheme.bed	
SWIFT_primer_bed	gs://theiagen-public-files/terra/theia	String Updated 2023-07-05 kk

Figure 18